



Mobile Dental Unit / Trolley CY-101


This sturdy cart features a large tray with rounded edges for holding instruments. It is equipped with a clean air, 600W, silent, oil free compressor with a stainless 12 liter tank.

*Easy drain.
Low maintenance.
No oil to change.*





Mobile Dental Unit / Trolley CY-101



Specifications

- 1 piece 3 way syringe water,air,mist w/600ml clean water bottle
- 3 pieces automatic handpiece selection controls
- 1 piece suction saliva absorber with 600ml collection bottle
- 1 piece low speed tubing 4 hole with LED control for LED polisher
- 1 piece high speed tubing 4 hole for LED high speed drill
- 1 piece LED ultrasonic scaler handpiece with 5 piece tips
- 1 piece LED high-speed push button handpiece
- 1 piece LED low speed polisher push button handpiece

distributed by



1.866.455.3050 | nuonimaging.ca

